

Product Datasheet

Anti-Herpes simplex virus, types 1&2 (TK) V0117

Product Name	Anti-Herpes simplex virus, types 1&2 (TK)
Catalog Number	V0117
Tested Applications	ELISA, WB, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 1 & 2
Specificity	Thymidine kinase
Shelf Life	Four years from date of manufacture