

Product Datasheet

Anti-Herpes simplex virus, type 2 (gG) V0125

Product Name	Anti-Herpes simplex virus, type 2 (gG)
Catalog Number	V0125
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 2. Non-reactive with HSV1, CMV, EBV, & VZV
Specificity	Glycoprotein G
Shelf Life	4 years