

# Product Datasheet

## Anti-Herpes simplex virus, types 1&2 (gD) V0194

Product Name	Anti-Herpes simplex virus, types 1&2 (gD)
Catalog Number	V0194
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 1 & 2
Specificity	Glycoprotein D
Shelf Life	Four years from date of manufacture