

Product Datasheet

Anti-Herpes simplex virus, types 1&2 (gD) V0195

Product Name	Anti-Herpes simplex virus, types 1&2 (gD)
Catalog Number	V0195
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 1 & 2
Specificity	Glycoprotein D
Matched Pairs	Pairs with V0197 in LF applications
Shelf Life	Four years from date of manufacture