

Product Datasheet

Anti-Herpes simplex virus, types 1&2 (gD) V0197

Product Name	Anti-Herpes simplex virus, types 1&2 (gD)
Catalog Number	V0197
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 1 & 2
Specificity	Glycoprotein D
Matched Pairs	Pairs with V0195 in LF applications
Shelf Life	4 years