

Product Datasheet

Anti-Epstein-Barr virus (Bam H1Z) V0256

Product Name	Anti-Epstein-Barr virus (Bam H1Z)
Catalog Number	V0256
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, IHC, ELISA, IP
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Early antigen (BAM HI Z) zeta protein
Shelf Life	4 years