

Product Datasheet

Anti-Epstein-Barr virus (EA-D) V0261

Product Name	Anti-Epstein-Barr virus (EA-D)
Catalog Number	V0261
Tested Applications	WB, IFA, IHC (weak)
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	“Early antigen-diffuse” (EA-D) of EBV
Shelf Life	Four years from date of manufacture