

Product Datasheet

Anti-Streptococcus pneumoniae V0491

Product Name	Anti-Streptococcus pneumoniae
Catalog Number	V0491
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Streptococcus pneumoniae common C-polysaccharide urinary antigen (backbone epitope)
Matched Pairs	ELISA and LF pairs (capture-conjugate): V0494-V0491, V0495-V0491
Shelf Life	4 years