

Product Datasheet

Anti-Rotavirus V0521

Product Name	Anti-Rotavirus
Catalog Number	V0521
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, LF, ELISA, ICC
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells. Additional cross reactivity data is available for this clone upon request.
Specificity	Recognizes the P41 capsid protein of monkey, porcine, & human isolates
Shelf Life	4 years