

Product Datasheet

Anti-Rotavirus V0521

| | |
|----------------------------|--|
| Product Name | Anti-Rotavirus |
| Catalog Number | V0521 |
| Description | Purified from mouse ascites by protein A chromatography |
| Species/Host | Mouse |
| Tested Applications | WB, ELISA, LF, ICC |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | +2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | ≥ 90 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells. Additional cross reactivity data is available for this clone upon request. |
| Specificity | Recognizes the P41 capsid protein of monkey, porcine, & human isolates |
| Shelf Life | 4 years |