

Product Datasheet

Anti-Rotavirus V0531

Product Name	Anti-Rotavirus
Catalog Number	V0531
Tested Applications	LF, ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with intact virus strains RRV, WA, SA11, & bovine. Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells
Specificity	Intact virus strains RRV, WA, SA11, & bovine
Matched Pairs	ELISA pairs (capture-conjugate): V0531-V0541, V0531-V0581
Shelf Life	Four years from date of manufacture