

Product Datasheet

Anti-Rotavirus V0541

Product Name	Anti-Rotavirus
Catalog Number	V0541
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with intact virus strains RRV, WA, SA11, & bovine. Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells
Specificity	Intact virus strains RRV, WA, SA11, & bovine
Matched Pairs	ELISA pairs (capture-probe): V0531-V0541
Shelf Life	4 years