

Product Datasheet

Anti-Respiratory Syncytial virus (fusion) V0651

Product Name	Anti-Respiratory Syncytial virus (fusion)
Catalog Number	V0651
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza A & B, Adenovirus, Para 1-3, Measles, & Mumps
Specificity	Fusion protein of RSV, types A & B
Matched Pairs	ELISA pairs (capture-conjugate) : V0631-V0651, V0651-V0635, V0635-V0651
Shelf Life	4 years