

Product Datasheet

Anti-Human CD4 Monoclonal Antibody, Azide Free Clone B-A1 854.030.005

Product Name Anti-Human CD4 Monoclonal Antibody, Azide Free Clone B-A1

Catalog Number 854.030.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications Functional assay, IHC, FC

Immunogen Human thymocytes and jurkat cell lines

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises a 55 kDa protein

Biological Activity Inhibits mixed lymphocyte culture

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells