

Product Datasheet

Anti-Cytomegalovirus (matrix 65 kDa) V0891

Product Name	Anti-Cytomegalovirus (matrix 65 kDa)
Catalog Number	V0891
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	65 kD late major matrix protein of CMV
Shelf Life	4 years