

Product Datasheet

Anti-Human CD4 Monoclonal Antibody, Unconjugated Clone B-A1 854.033.020

Product Name	Anti-Human CD4 Monoclonal Antibody, Unconjugated Clone B-A1
Catalog Number	854.033.020
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IHC, FC
Immunogen	Human thymocytes and jurkat cell lines
Brand	Diaclone
Form/Appearance	Liquid
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognises a 55 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells