

# Product Datasheet

## Anti-E. coli V1011

|                     |   |
|---------------------|---|
| Product Name        | Anti-E. coli  |
| Catalog Number      | V1011   |
| Tested Applications | ELISA assay   |
| Brand               | ViroStat  |
| Form/Appearance     | Purified from mouse ascites by protein A chromatography   |
| Concentration       | ≥ 1 mg/ml   |
| Storage             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.                                       |
| Isotype             | IgG1  |
| Clonality           | Monoclonal  |
| Purity              | ≥ 90%   |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide   |
| Cross Reactivity    | Reactive with at least serotypes O18, O44, O112, & O125.<br>Non-reactive with other members of Enterobacteriaceae |
| Specificity         | E. coli serotypes   |
| Shelf Life          | Four years from date of manufacture   |