

# Product Datasheet

## Anti-E. coli V1011

Product Name	Anti-E. coli
Catalog Number	V1011
Tested Applications	ELISA assay
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with at least serotypes O18, O44, O112, & O125. Non-reactive with other members of Enterobacteriaceae
Specificity	E. coli serotypes
Shelf Life	Four years from date of manufacture