

Product Datasheet

Anti-E. coli O157 V1023

Product Name	Anti-E. coli O157
Catalog Number	V1023
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	“O” antigen of E. coli serotype O157
Matched Pairs	ELISA pairs (capture-conjugate) : V1023-V1022, V1026-V1023
Shelf Life	Four years from date of manufacture