

Product Datasheet

Anti-Influenza B virus (NP) V1131

Product Name	Anti-Influenza B virus (NP)
Catalog Number	V1131
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Influenza B nucleoprotein
Matched Pairs	LF pairs (capture-conjugate) : V1161-V1131, V1181-V1131, V1171-V1131
Shelf Life	4 years