

# Product Datasheet

## Anti-Parainfluenza virus type 1 V1226

Product Name	Anti-Parainfluenza virus type 1
Catalog Number	V1226
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Fusion protein of Parainfluenza virus type 1
Matched Pairs	Pairs in ELISA with V1224
Shelf Life	Four years from date of manufacture