

Product Datasheet

Anti-Parainfluenza virus type 1 V1231

Product Name	Anti-Parainfluenza virus type 1
Catalog Number	V1231
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with types 2 & 3 by indirect IFA. Reactive with Sendai virus
Specificity	Parainfluenza virus type 1
Shelf Life	Four years from date of manufacture