

Product Datasheet

Anti-Influenza A virus (NP) V1341

Product Name	Anti-Influenza A virus (NP)
Catalog Number	V1341
Tested Applications	ELISA, IFA, Cyto
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with all A strains
Specificity	Influenza A – specific for nucleoprotein
Shelf Life	Four years from date of manufacture