

## Product Datasheet

### Anti-Adenovirus V1431

<b>Product Name</b>	Anti-Adenovirus
<b>Catalog Number</b>	V1431
<b>Description</b>	Purified from mouse ascites by protein A chromatography
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with 34 serotypes of Adenovirus including types 40 & 41. Non-reactive with Influenza A, Influenza B, RSV, Para 1, 2, & 3, Mycoplasma pneumoniae, H. pylori and mammalian cells. Additional cross reactivity data is available upon request.
<b>Specificity</b>	Adenovirus hexon group antigen
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate) : V1431-V1411, V1431-V1461, V1431-V1481, V1411- V1431, V1421-V1461, V1421-V1411
<b>Shelf Life</b>	4 years