

Product Datasheet

Anti-Adenovirus V1461

Product Name	Anti-Adenovirus
Catalog Number	V1461
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with at least 21 serotypes of Adenovirus including types 40 & 41. Non-reactive with Influenza A, Influenza B, RSV, Para 1, 2, & 3, Mycoplasma pneumoniae, H. pylori and mammalian cells
Specificity	Adenovirus hexon group antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V1431-V1411, V1431-V1461, V1431-V1481, V1411- V1431, V1421-V1461, V1421-V1411
Shelf Life	4 years