

## Product Datasheet

### Anti-Chlamydia (LPS) V1646

<b>Product Name</b>	Anti-Chlamydia (LPS)
<b>Catalog Number</b>	V1646
<b>Description</b>	Purified from mouse ascites by protein A chromatography
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with <i>C. trachomatis</i> , <i>C. pneumoniae</i> , <i>C. psittaci</i>
<b>Specificity</b>	LPS
<b>Matched Pairs</b>	LPS ELISA pairs (capture-conjugate): V1646-V1644
<b>Shelf Life</b>	4 years