

Product Datasheet

Anti-Hepatitis B virus (HBsAg) V1836

Product Name	Anti-Hepatitis B virus (HBsAg)
Catalog Number	V1836
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, IHC, ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with normal human serum
Specificity	Surface antigen; subtypes ad & ay (clone developed using native antigen)
Matched Pairs	ELISA pairs (capture-conjugate) : V1836-V1837, V1836-V1814
Shelf Life	4 years