

Product Datasheet

Anti-Hepatitis C virus (core) V1851

Product Name	Anti-Hepatitis C virus (core)
Catalog Number	V1851
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with genotypes 1a, 1b, and 2a
Specificity	Core, 29-39 aa
Shelf Life	4 years