

Product Datasheet

Anti-Human CD7 Monoclonal Antibody, Azide Free Clone B-F12 854.060.005

Product Name	Anti-Human CD7 Monoclonal Antibody, Azide Free Clone B-F12
Catalog Number	854.060.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Chronic Lymphoid Leukemia cells (CLL)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises CD7, a 40 kD type I transmembrane glycoprotein expressed on peripheral blood T lymphocytes, NK-cells and hematopoietic progenitors.
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells