

Product Datasheet

Anti-Hepatitis B virus (HBeAg) V1892

Product Name	Anti-Hepatitis B virus (HBeAg)
Catalog Number	V1892
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Hepatitis B “e” antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V1892-V1893, V1892-V1894
Shelf Life	4 years