

Product Datasheet

Anti-Treponema pallidum V2121

Product Name	Anti-Treponema pallidum
Catalog Number	V2121
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with a 30-32 kDa of T. pallidum. Cross reactive with T. reiter. Non-reactive with Borrelia burgdorferi
Specificity	30-32 kDa of T. pallidum
Shelf Life	Four years from date of manufacture