

# Product Datasheet

## Anti-Treponema pallidum V2124

Product Name	Anti-Treponema pallidum
Catalog Number	V2124
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Borrelia burgdorferi
Specificity	TmpA of T. pallidum
Shelf Life	Four years from date of manufacture