

Product Datasheet

Anti-Treponema pallidum V2126

Product Name	Anti-Treponema pallidum
Catalog Number	V2126
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Borrelia burgdorferi
Specificity	T. pallidum p47
Matched Pairs	V2126 as capture + V2125 as conjugate
Shelf Life	4 years