

Product Datasheet

Anti-Treponema pallidum V2128

Product Name	Anti-Treponema pallidum
Catalog Number	V2128
Tested Applications	ELISA, WB, IFA (does not stain the entire bacterial cell)
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Borrelia burgdorferi
Specificity	T. pallidum p47
Shelf Life	Four years from date of manufacture