

Product Datasheet

Anti-Ebola virus (NP) V2601

| | |
|----------------------------|---|
| Product Name | Anti-Ebola virus (NP) |
| Catalog Number | V2601 |
| Tested Applications | LF, ELISA, WB |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Reactive with the NP of Sudan & Zaire sequences & VLPs in ELISA & WB |
| Specificity | NP of Sudan & Zaire sequences |
| Matched Pairs | LF pairs (capture-conjugate) : V2601-V2602 |
| Shelf Life | Four years from date of manufacture |