

Product Datasheet

Anti-Ebola virus (NP) V2601

Product Name	Anti-Ebola virus (NP)
Catalog Number	V2601
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, LF
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the NP of Sudan & Zaire sequences & VLPs in ELISA & WB
Specificity	NP of Sudan & Zaire sequences
Matched Pairs	LF pairs (capture-conjugate) : V2601-V2602
Shelf Life	4 years