

# Product Datasheet

## Anti-Yellow Fever virus (NS1) V3573

Product Name	Anti-Yellow Fever virus (NS1)
Catalog Number	V3573
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive in IFA vs TBE virus, WNV & JEV
Specificity	NS1
Shelf Life	Four years from date of manufacture