

## Product Datasheet

### Anti-Chikungunya virus V3585

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Chikungunya virus   |
| <b>Catalog Number</b>      | V3585  |
| <b>Tested Applications</b> | ELISA, IFA (membrane staining)   |
| <b>Brand</b>               | ViroStat   |
| <b>Form/Appearance</b>     | Purified from mouse ascites by protein A chromatography                            |
| <b>Concentration</b>       | ≥ 1 mg/ml  |
| <b>Storage</b>             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.        |
| <b>Isotype</b>             | IgG2b  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | ≥ 90%  |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide              |
| <b>Cross Reactivity</b>    | Non-reactive with other Alphaviruses including EEE, WEE, VEE and other Arboviruses |
| <b>Specificity</b>         | Chikungunya virus  |
| <b>Shelf Life</b>          | Four years from date of manufacture  |