

# Product Datasheet

## Anti-West Nile Virus (Env) V3641

Product Name	Anti-West Nile Virus (Env)
Catalog Number	V3641
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	A table of cross reactivity with flaviviruses is available upon request
Specificity	Envelope protein
Shelf Life	Four years from date of manufacture