

Product Datasheet

Anti-Enterovirus 71 V3657

Product Name	Anti-Enterovirus 71
Catalog Number	V3657
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Coxsackie A9, Coxsackie B2, ECHO 30, EV70 and Polio, type 2
Specificity	VP1 of Enterovirus type 71
Shelf Life	Four years from date of manufacture