

Product Datasheet

Anti-West Nile Virus (Env) V3671

| | |
|----------------------------|--|
| Product Name | Anti-West Nile Virus (Env) |
| Catalog Number | V3671 |
| Description | Purified from mouse ascites by protein A chromatography |
| Species/Host | Mouse |
| Tested Applications | ELISA, IFA |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | +2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | ≥ 90 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | A table of cross reactivity with flaviviruses is available upon request |
| Specificity | Envelope protein |
| Shelf Life | 4 years |