

Product Datasheet

Anti-Human CD10 Monoclonal Antibody, Azide Free Clone B-E3 854.090.005

Product Name	Anti-Human CD10 Monoclonal Antibody, Azide Free Clone B-E3
Catalog Number	854.090.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay, IHC, FC
Immunogen	Daudi and Raji cell lines
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Common Acute Lymphoblastic Leukemia Antigen (CALLA), a 100 kDa protein
Biological Activity	Cytotoxic In association with rabbit complement on CD10+ cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells