

Product Datasheet

Anti-Powassan Virus V3683

| | |
|---------------------|---|
| Product Name | Anti-Powassan Virus |
| Catalog Number | V3683 |
| Tested Applications | ELISA, IFA |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with West Nile virus, Eastern Equine Encephalitis virus, Western Equine Encephalitis virus, St. Louis Encephalitis virus and La Crosse virus |
| Specificity | Powassan virus envelope III |
| Shelf Life | Four years from date of manufacture |