

Product Datasheet

Anti-Human CD10 Monoclonal Antibody, FITC Conjugated Clone MEM-78 854.091.010

Product Name	Anti-Human CD10 Monoclonal Antibody, FITC Conjugated Clone MEM-78
Catalog Number	854.091.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Pre-B cell line NALM-6
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA)
Storage	Store at +2-8°C protected from light. DO NOT FREEZE.
Note	For flow cytometry we recommend using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension
Isotype	IgG1
Clonality	Monoclonal
Source	Pre-B cell line NALM-6
Specificity	Reacts with CD10 antigen (CALLA - Common acute lymphatic leukemia antigen), a 100 kDa type II integral membrane protein expressed on uncommitted lymphoid precursors