

## Product Datasheet

### Anti-Coronavirus (NP) V3879

**Product Name** Anti-Coronavirus (NP)

**Catalog Number** V3879

**Tested Applications** ELISA, IFA

**Brand** ViroStat

**Form/Appearance** Purified from mouse ascites by protein A chromatography

**Concentration**  $\geq 1$  mg/ml

**Storage** 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

**Isotype** IgG1

**Clonality** Monoclonal

**Purity**  $\geq 90\%$

**Buffer** 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

**Cross Reactivity** Reactive with NP of native SARS, MERS by IFA, & COVID-19 by ELISA. Non-reactive with Coronavirus strains 229E, NL63, & OC43

**Specificity** NP of native SARS, MERS, & COVID-19

**Shelf Life** Four years from date of manufacture