

Product Datasheet

Anti-Coronavirus (NP) V3879

Product Name	Anti-Coronavirus (NP)
Catalog Number	V3879
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with NP of native SARS, MERS by IFA, & COVID-19 by ELISA. Non-reactive with Coronavirus strains 229E, NL63, & OC43
Specificity	NP of native SARS, MERS, & COVID-19
Shelf Life	4 years