

## Product Datasheet

### Anti-Listeria monocytogenes V4211

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Listeria monocytogenes  |
| <b>Catalog Number</b>      | V4211  |
| <b>Description</b>         | Purified from mouse ascites by protein A chromatography                        |
| <b>Species/Host</b>        | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Brand</b>               | ViroStat   |
| <b>Form/Appearance</b>     | Purified from mouse ascites by protein A chromatography                        |
| <b>Concentration</b>       | ≥ 1 mg/ml  |
| <b>Storage</b>             | +2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | ≥ 90 %   |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide          |
| <b>Cross Reactivity</b>    | Reactive with a membrane antigen of L. monocytogenes                           |
| <b>Specificity</b>         | Membrane antigen of L. monocytogenes   |
| <b>Matched Pairs</b>       | ELISA pairs (capture-conjugate): V4211-V4212                                   |
| <b>Shelf Life</b>          | 4 years  |