

Product Datasheet

Anti-Metapneumovirus V4831

Product Name	Anti-Metapneumovirus
Catalog Number	V4831
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	IHC, IFA
Alternative Names	HMPV
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Respiratory Syncytial Virus, Influenza A, Influenza B, Adenovirus, & Parainfluenza viruses 1, 2, & 3
Specificity	Matrix protein of human metapneumovirus
Shelf Life	4 years