

# Product Datasheet

## Anti-Metapneumovirus V4841

Product Name	Anti-Metapneumovirus
Catalog Number	V4841
Tested Applications	ELISA, IFA
Alternative Names	HMPV
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Respiratory Syncytial Virus, Influenza A, Influenza B, Adenovirus, & Parainfluenza viruses 1, 2, & 3
Specificity	Matrix protein of human metapneumovirus
Shelf Life	Four years from date of manufacture