

Product Datasheet

Anti-Norovirus (GI) V4922

| | |
|----------------------------|--|
| Product Name | Anti-Norovirus (GI) |
| Catalog Number | V4922 |
| Description | Purified from mouse ascites by protein A chromatography |
| Species/Host | Mouse |
| Tested Applications | WB, ELISA, LF |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | +2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with GII |
| Specificity | GI capsid protein |
| Matched Pairs | ELISA pairs (capture-conjugate) : V4922-V4924. LF IA pairs (capture-conjugate) : V4925-V4922 |
| Shelf Life | 4 years |